



RevoFit[®]

Design Fundamentals


How to choose the right design based on your patient indications.

There are 3 proven design options to successfully build RevoFit[®] adjustability in just about any device:

- **Panel** design
- **Gap** design
- **Hinge** design

Panel

**BEST FOR TARGETED
COMPRESSION / EXPANSION**



Solves For

- ✓ Bone levering
- ✓ Rotation control
- ✓ Bone or muscle lock
- ✓ Volume fluctuations
- ✓ Donn/doffing issues

- **Ideal panel locations:**

Pressure tolerant areas

- **Panel size:**

Customize shape and size to users anatomy, similar to the pads you would install in a rigid device.


- **Pro-tip:**

Test panel location assumptions in the diagnostic phase.



Gap

**BEST FOR GLOBAL
COMPRESSION / EXPANSION**



Solves For

- ✓ Bone or muscle lock
- ✓ Volume fluctuations
- ✓ Donn/doffing issues
- ✓ Device instability
- ✓ Expansion needs large distal ends

• Ideal panel locations:


- For TF: Across proximal posterior opening
- For TT: Posterior medial/lateral corners

• Pro-tip:

Gaps only close at the proximal level. The longer the gap extends distally, the easier the proximal portion will close.

Hinge

**BEST FOR
DONN/DOFF AND SUSPENSION**



Solves For

- ✓ Volume fluctuations
- ✓ Donn/doffing issues
- ✓ Device instability
- ✓ Primary suspension
- ✓ Expansion needs for large distal ends

• Ideal panel locations:

- Posterior aspect for TT
- Medial aspect for supracondylar suspension

• Safety Tip:

Creating a posterior proximal brim maintains structural integrity and ensures safety if the hinge were to open unexpectedly.