Instructions For Use

Meracus®

Art. No. 29000

Instructions For Use

EN	Meracus® Prosthetic Foot	Page 4
SV	Meracus® Protesfot	Page 8
NO	Meracus® Protese fot	Page 12
DA	Meracus® Fodprotese	Page 16
FI	Meracus® Jalkateräproteesi	Page 20
ES	Meracus® Pie Protésico	Page 24
DE	Meracus® Prothesenfuß	Page 28
NL	Meracus® Prothesevoet	Page 32
FR	Meracus® Pied prothétique	Page 36
IT	Meracus® Piede protesico	Page 40

EN Meracus Prosthetic Foot

Product Description

Meracus is an energy-storing, dynamic response prosthetic foot that consists of a pyramid adapter, J-shaped composite keel and a full-length composite heel plate.

Meracus has dynamic properties that simulate the anatomical foot by providing shock absorption during heel strike, smooth transition to a stable stance phase, smooth roll-off of the forefoot and efficient energy return during push-off.

Intended Use

The prosthetic foot is intended to replace the function of the amputated foot and ankle. The prosthetic foot is suitable for moderate to active use, e.g., walking, hiking and occasional running. Single patient - multiple use.

Indications

Lower extremity limb loss or limb difference.

Contraindications

This device may not be suitable for very low active individuals or those engaged in competitive sports or other very high impact activities, nor for patients not able to read and understand or follow user instructions and guidelines.

Limitations

Weight limitations: for the larger sizes 150kg (330lb) and for the smallest size 100kg (220lb). See detailed information under Product Selection.

Patient assessment

Adult individuals that has the ability or potential for ambulation with variable cadence. Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

Material Specification

- Carbon fiber and glass fiber in an epoxy matrix
- Titanium
- Stainless Steel
- EPDM



Functional life

The prosthetic foot has undergone a dynamic test in two million cycles, in accordance with standard ISO 22675. The expected functional life of the product is approximately 2-3 years, depending on the patient's activity level and adherence to guidelines.



Product Selection Guide

Weight, kg	≤60	≤80	≤100	≤125	≤150
Weight, lbs	≤132	≤176	≤220	≤275	≤330
	Standard Pyramid Adapter			Extended Pyramid Adapter	
Size 23	23P3	23P4	23P5	Not available	
Build height	125mm	125mm	125mm	Not available	
Size 24	24P3	24P4	24P5	24P6	Not available
Build height	125mm	125mm	125mm	133mm	Not available
Size 25	25P3	25P4	25P5	25P6	25P7
Build height	132mm	132mm	132mm	140mm	140mm
Size 26	26P3	26P4	26P5	26P6	26P7
Build height	132mm	132mm	132mm	140mm	140mm
Size 27	27P3	27P4	27P5	27P6	27P7
Build height	136mm	136mm	136mm	144mm	144mm
Size 28	Not available	28P4	28P5	28P6	28P7
Build height	Not available	136mm	136mm	144mm	144mm

Our selection chart is based on testing at each weight level. If a stiffer foot is preferred, you may move up one weight category. When choosing a product, also consider the user's daily activities, such as heavy lifting in connection with work or leisure activities/sports. Please note that changing from P5 to P6 results in an 8mm increase in build height.

Example:

EN Meracus Prosthetic Foot

Initial bench alignment

The prosthetic foot must be aligned by a prosthetist or under supervision of a prosthetist, or equivalent medical professional. Standard bench adjustment techniques can be applied for Meracus and are only an initial guide. Place the prosthetic foot on a flat surface with a 10mm block under the heel to simulate the heel height. The plumb line shall pass through the center of the prosthetic socket and to a point at approximately 1/3 of the whole length of the foot (fig. 1).

Use a plumbline or laser to divide the socket into equal portions Medial and Lateral. The plumbline should go through the center of the pyramid attachment and bisect to heel (fig. 2).


Fig. 1


Fig. 2

Static and dynamic alignment

Meracus is designed with a heel height of 10mm. The static and dynamic alignment of the foot should be done with the shoes (or shoe type) the patient wears most often. Shoes with softer or firmer soles will affect the functional alignment and performance of the foot. The static adjustment is done in a standing position while the patient puts weight on the prosthesis to make all the necessary corrections in height, M-L and A-P alignment.

The dynamic alignment is done as the patient walks, to fine-tune the transition from initial heel contact to energy return, and toe-off. If the progression of the lower leg is too fast from heel to midstance, consider shifting the socket posteriorly and/ or plantarflexing the foot slightly. If the heel is too soft, the opposite situation may be considered; an anterior shift of the socket and/ or dorsal flexion of the foot, to correct a slowed progression of the lower leg and heel-to-toe movement.

The same principles apply to the keel if lower leg progression is delayed from midstance to toe-off; consider an anterior shift of the socket anterior and/ or dorsiflexing the foot slightly. If the progression is too quick from midstance to pre-swing, consider shifting the socket posterior and/ or plantarflexing the foot slightly.

Once the alignment process is complete all screws should have thread locking compound (Loctite) applied and be torqued to the manufacturers recommended specifications. Once the alignment process is complete all screws should have thread locking compound (Loctite) applied and be torqued to the manufacturers recommended specifications.

Information

Any serious incident that may occur in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and / or patient is established.

Disposal

To be disposed in accordance with local legislations.

Warnings

The prosthetist or equivalent medical professional must review the following information directly with the user:

- Meracus is tested and approved according to ISO 22675 and should be used together with standardized adult-size components. It is the fitter's responsibility to ensure safe combinations of componentry in the final prosthetic device.
- The prosthetic foot is water resistant, but do not let debris or liquid stay inside the foot shell. If the prosthesis encounters either, (such as sand, salt water or chlorinated pools) remove the foot shell by proper usage of a foot shell removal tool to protect your fingers. Take off the protective Spectra sock, rinse it with clean water, and let it dry completely. Inspect the foot to make sure no debris remains.
- The prosthetic foot should be removed from the foot shell and inspected by clinicians every six months to detect any signs of abnormal wear and to check that the screws are secured.
- If any changes or modifications are made to the product that are not described in the Professional Instructions, the product responsibility passes entirely to the person making the changes.
- The product should not be used for activities where there is a risk of being exposed to exceptional strain and/or weight. Not suitable for high impact sports.
- Evaluate the ability to drive or operate heavy machinery. It is important to consider that a prosthetic device may create hazards not present with an intact limb.
- The prosthetic foot must always be used with a Spectra sock and a foot shell. The Spectra sock protects the foot shell and composite from wear.
- The prosthetic foot, when used with only the foot shell and no shoe, can be slippery—particularly in wet environments such as during showering or bathing.
- If the properties of the prosthetic foot change, or if an abnormal noise is heard from the prosthetic foot, the fitter must be contacted immediately, and the foot must not be used.
- Use the pyramid protector during the designated trial procedures. Ensure that the protector is removed before the patient leaves the clinic or if the foot is intended to be taken off-site.

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Camp Scandinavia AB