

## PROSTHETIC SHEATHS

Sheaths control friction, adjust thickness, provide strength and aid skin care. 100% stretch nylon sheaths come in dusk and white colors, four widths which correspond to prosthetic sock sizes. Super Stretch Sheaths have Lycra added to the nylon to accommodate an even greater range and more conforming fit, enhance flexion, movement and improve suspension. Knit-Rite Sheaths have transfers on both sheaths of the pair to help patients position the sheath to maximize comfort.

The Silver Sheath with X-Static® is anti-odor. It also wicks moisture away from the skin and aids in thermal consistency.



PROPER CARE: Sheaths may be hand or machine washed. Use soap or mild detergent in warm water. NO BLEACH. Rinse thoroughly. Air dry overnight or tumble dry low.

SILVER	DUSK	WHITE	SOCK SIZE	LENGTH
<b>SILVER SHEATH — 91% NYLON / 9% PURE SILVER X-STATIC®. SOLD PER BOX (3 PAIR).</b>				
<b>2H-PXNA-XS</b>	—	—	Narrow X-Short #B, #O, #1	8"-10"
<b>2H-PXNA-SH</b>	—	—	Narrow Short #B, #O, #1	10"-14"
<b>2H-PXNA-MD</b>	—	—	Narrow Medium #B, #O, #1	16"-20"
<b>2H-PXRG-XS</b>	—	—	Regular X-Short #1, #2	8"-10"
<b>2H-PXRG-SH</b>	—	—	Regular Short #1, #2	10"-14"
<b>2H-PXRG-MD</b>	—	—	Regular Medium #1, #2	16"-20"
<b>2H-PXWD-SH</b>	—	—	Wide Short #3 and above	10"-14"

**SUPER STRETCH SHEATH — 98% NYLON, 2% LYCRA. SOLD PER BOX (3 PAIR).**

—	<b>2H-PSBR-SH</b>	<b>2H-PSWR-SH</b>	Regular Short #1, #2	10"-14"+
—	<b>2H-PSBR-MD</b>	<b>2H-PSWR-MD</b>	Regular Medium #1, #2	16"-20"+

**PROSTHETIC SHEATH — 100% STRETCH NYLON, SIZED TO SOCKS. SOLD PER BOX (3 PAIR).**

—	<b>2H-PBCH-SH</b>	<b>2H-PWCH-SH</b>	Child Short #A, Lycra spandex top	6"-8"
—	<b>2H-PBCH-12</b>	<b>2H-PWCH-12</b>	Child 12" #A, non-elastic top	12"
—	<b>2H-PBNA-SH</b>	<b>2H-PWNA-SH</b>	Narrow Short #B, #O, #1	10"-14"
—	<b>2H-PBNA-MD</b>	<b>2H-PWNA-MD</b>	Narrow Medium #B, #O, #1	16"-20"
—	<b>2H-PBNA-LG</b>	<b>2H-PWNA-LG</b>	Narrow Long #B, #O, #1	22"-28"
—	<b>2H-PBRG-XS</b>	<b>2H-PWRG-XS</b>	Regular X-Short #1, #2	8"-10"
—	<b>2H-PBRG-SH</b>	<b>2H-PWRG-SH</b>	Regular Short #1, #2	10"-14"
—	<b>2H-PBRG-MD</b>	<b>2H-PWRG-MD</b>	Regular Medium #1, #2	16"-20"
—	<b>2H-PBRG-LG</b>	<b>2H-PWRG-LG</b>	Regular Long #1, #2	22"-28"
—	<b>2H-PBWD-SH</b>	<b>2H-PWWD-SH</b>	Wide Short #3 and above	10"-14"
—	<b>2H-PBWD-MD</b>	<b>2H-PWWD-MD</b>	Wide Medium #3 and above	16"-20"
—	<b>2H-PBWD-LG</b>	<b>2H-PWWD-LG</b>	Wide Long #3 and above	22"-28"

**HOLE-IN-TOE SHEATH — 100% STRETCH NYLON. SOLD PER PAIR.**

—	<b>2H-PBRH-SH</b>	<b>2H-PWRH-SH</b>	Regular Short #1, #2	10"-14"
—	<b>2H-PBWH-SH</b>	<b>2H-PWWH-SH</b>	Wide Short #3 and above	10"-14"

