

# RUSH HIPRO

Kategorivelger i forhold til aktivitetsnivå, brukervekt og størrelse																		
Vekt	0 - 48 kg		49 - 64 kg		64 - 79 kg		80 - 95 kg		96 - 111 kg		112 - 127 kg		128 - 143 kg		144 - 159 kg		160 - 166 kg	
	Alle	High Impact	Alle	High Impact	Alle	High Impact	Alle	High Impact	Alle	High Impact	Alle	High Impact	Alle	High Impact	Alle	High Impact	Alle	High Impact
22 cm	1	2	2	3	3	4	4	5	5	6	6	-	-	-	-	-	-	-
23 cm	1	2	2	3	3	4	4	5	5	6	6	-	-	-	-	-	-	-
24 cm	-	-	2	3	3	4	4	5	5	6	6	-	-	-	-	-	-	-
25 cm	-	-	2	3	3	4	4	5	5	6	6	7	7	-	-	-	-	-
26 cm	-	-	2	3	3	4	4	5	5	6	6	7	7	8	8	-	-	-
27 cm	-	-	-	-	3	4	4	5	5	6	6	7	7	8	8	9	9	-
28 cm	-	-	-	-	3	4	4	5	5	6	6	7	7	8	8	9	9	-
29 cm	-	-	-	-	3	4	4	5	5	6	6	7	7	8	8	9	9	-

Art.nr.	Modell	Lev.art.nr.
759022100	RUSH HIPRO 22cm cat 1	HIP-22-1-RU
759022200	RUSH HIPRO 22cm cat 2	HIP-22-2-RU
759022300	RUSH HIPRO 22cm cat 3	HIP-22-3-RU
759022400	RUSH HIPRO 22cm cat 4	HIP-22-4-RU
759022500	RUSH HIPRO 22cm cat 5	HIP-22-5-RU
759023100	RUSH HIPRO 23cm cat 1	HIP-23-1-RU
759023200	RUSH HIPRO 23cm cat 2	HIP-23-2-RU
759023300	RUSH HIPRO 23cm cat 3	HIP-23-3-RU
759023400	RUSH HIPRO 23cm cat 4	HIP-23-4-RU
759023500	RUSH HIPRO 23cm cat 5	HIP-23-5-RU
759024200	RUSH HIPRO 24cm cat 2	HIP-24-2-RU
759024300	RUSH HIPRO 24cm cat 3	HIP-24-3-RU
759024400	RUSH HIPRO 24cm cat 4	HIP-24-4-RU
759024500	RUSH HIPRO 24cm cat 5	HIP-24-5-RU
759024600	RUSH HIPRO 24cm cat 6	HIP-24-6-RU
759025200	RUSH HIPRO 25cm cat 2	HIP-25-2-RU
759025300	RUSH HIPRO 25cm cat 3	HIP-25-3-RU
759025400	RUSH HIPRO 25cm cat 4	HIP-25-4-RU
759025500	RUSH HIPRO 25cm cat 5	HIP-25-5-RU
759025600	RUSH HIPRO 25cm cat 6	HIP-25-6-RU
759025700	RUSH HIPRO 25cm cat 7	HIP-25-7-RU
759026200	RUSH HIPRO 26cm cat 2	HIP-26-2-RU
759026300	RUSH HIPRO 26cm cat 3	HIP-26-3-RU
759026400	RUSH HIPRO 26cm cat 4	HIP-26-4-RU
759026500	RUSH HIPRO 26cm cat 5	HIP-26-5-RU
759026600	RUSH HIPRO 26cm cat 6	HIP-26-6-RU
759026700	RUSH HIPRO 26cm cat 7	HIP-26-7-RU
759027300	RUSH HIPRO 27cm cat 3	HIP-27-3-RU
759027400	RUSH HIPRO 27cm cat 4	HIP-27-4-RU
759027500	RUSH HIPRO 27cm cat 5	HIP-27-5-RU
759027600	RUSH HIPRO 27cm cat 6	HIP-27-6-RU
759027700	RUSH HIPRO 27cm cat 7	HIP-27-7-RU
759027800	RUSH HIPRO 27cm cat 8	HIP-27-8-RU
759027900	RUSH HIPRO 27cm cat 9	HIP-27-9-RU
759028300	RUSH HIPRO 28cm cat 3	HIP-28-3-RU
759028400	RUSH HIPRO 28cm cat 4	HIP-28-4-RU
759028500	RUSH HIPRO 28cm cat 5	HIP-28-5-RU
759028600	RUSH HIPRO 28cm cat 6	HIP-28-6-RU
759028700	RUSH HIPRO 28cm cat 7	HIP-28-7-RU
759028800	RUSH HIPRO 28cm cat 8	HIP-28-8-RU
759028900	RUSH HIPRO 28cm cat 9	HIP-28-9-RU
759029300	RUSH HIPRO 29cm cat 3	HIP-29-3-RU
759029400	RUSH HIPRO 29cm cat 4	HIP-29-4-RU
759029500	RUSH HIPRO 29cm cat 5	HIP-29-5-RU
759029600	RUSH HIPRO 29cm cat 6	HIP-29-6-RU
759029700	RUSH HIPRO 29cm cat 7	HIP-29-7-RU
759029800	RUSH HIPRO 29cm cat 8	HIP-29-8-RU
759029900	RUSH HIPRO 29cm cat 9	HIP-29-9-RU