

RUSH ROGUE 2

Kategoriøvelser i forhold til aktivitetsnivå, brukervekt og størrelse																					
Vekt	0-48 kg		49-64 kg		64-79 kg		80-95 kg		96-111 kg		112-127 kg		128-143 kg		144-159 kg		160-166 kg				
	Alle	High Impact	Alle	High Impact	Alle	High Impact	Alle	High Impact	Alle	High Impact	Alle	High Impact	Alle	High Impact	Alle	High Impact	Alle	High Impact			
22 cm	1	2	2	3	3	4	4	5	5	6	6	-	-	-	-	-	-	-	-		
23 cm	1	2	2	3	3	4	4	5	5	6	6	-	-	-	-	-	-	-	-		
24 cm	-	-	2	3	3	4	4	5	5	6	6	-	-	-	-	-	-	-	-		
25 cm	-	-	2	3	3	4	4	5	5	6	6	7	7	-	-	-	-	-	-		
26 cm	-	-	2	3	3	4	4	5	5	6	6	7	7	8	8	-	-	-	-		
27 cm	-	-	-	-	3	4	4	5	5	6	6	7	7	8	8	9	9	-	-		
28 cm	-	-	-	-	3	4	4	5	5	6	6	7	7	8	8	9	9	-	-		
29 cm	-	-	-	-	3	4	4	5	5	6	6	7	7	8	8	9	9	-	-		

Art.nr.	Modell	Lesart.nr.
759322100	RUSH ROGUE 2 22cm cat 1	ROG2-22-1-RU
759322110	RUSH ROGUE 2 22cm Cat 1 SL	ROG2-22-1-SL
759322120	RUSH ROGUE 2 22cm Cat 1 SR	ROG2-22-1-SR
759322200	RUSH ROGUE 2 22cm cat 2	ROG2-22-2-RU
759322210	RUSH ROGUE 2 22cm Cat 2 SL	ROG2-22-2-SL
759322220	RUSH ROGUE 2 22cm Cat 2 SR	ROG2-22-2-SR
759322300	RUSH ROGUE 2 22cm cat 3	ROG2-22-3-RU
759322310	RUSH ROGUE 2 22cm Cat 3 SL	ROG2-22-3-SL
759322320	RUSH ROGUE 2 22cm Cat 3 SR	ROG2-22-3-SR
759322400	RUSH ROGUE 2 22cm cat 4	ROG2-22-4-RU
759322410	RUSH ROGUE 2 22cm Cat 4 SL	ROG2-22-4-SL
759322420	RUSH ROGUE 2 22cm Cat 4 SR	ROG2-22-4-SR
759322500	RUSH ROGUE 2 22cm cat 5	ROG2-22-5-RU
759322510	RUSH ROGUE 2 22cm Cat 5 SL	ROG2-22-5-SL
759322520	RUSH ROGUE 2 22cm Cat 5 SR	ROG2-22-5-SR
759323100	RUSH ROGUE 2 23cm cat 1	ROG2-23-1-RU
759323110	RUSH ROGUE 2 23cm Cat 1 SL	ROG2-23-1-SL
759323120	RUSH ROGUE 2 23cm Cat 1 SR	ROG2-23-1-SR
759323200	RUSH ROGUE 2 23cm cat 2	ROG2-23-2-RU
759323210	RUSH ROGUE 2 23cm Cat 2 SL	ROG2-23-2-SL
759323220	RUSH ROGUE 2 23cm Cat 2 SR	ROG2-23-2-SR
759323300	RUSH ROGUE 2 23cm cat 3	ROG2-23-3-RU
759323310	RUSH ROGUE 2 23cm Cat 3 SL	ROG2-23-3-SL
759323320	RUSH ROGUE 2 23cm Cat 3 SR	ROG2-23-3-SR
759323400	RUSH ROGUE 2 23cm cat 4	ROG2-23-4-RU
759323410	RUSH ROGUE 2 23cm Cat 4 SL	ROG2-23-4-SL
759323420	RUSH ROGUE 2 23cm Cat 4 SR	ROG2-23-4-SR
759323500	RUSH ROGUE 2 23cm cat 5	ROG2-23-5-RU
759323510	RUSH ROGUE 2 23cm Cat 5 SL	ROG2-23-5-SL
759323520	RUSH ROGUE 2 23cm Cat 5 SR	ROG2-23-5-SR
759324200	RUSH ROGUE 2 24cm cat 2	ROG2-24-2-RU
759324210	RUSH ROGUE 2 24cm Cat 2 SL	ROG2-24-2-SL
759324220	RUSH ROGUE 2 24cm Cat 2 SR	ROG2-24-2-SR
759324300	RUSH ROGUE 2 24cm cat 3	ROG2-24-3-RU
759324310	RUSH ROGUE 2 24cm Cat 3 SL	ROG2-24-3-SL
759324320	RUSH ROGUE 2 24cm Cat 3 SR	ROG2-24-3-SR
759324400	RUSH ROGUE 2 24cm cat 4	ROG2-24-4-RU
759324410	RUSH ROGUE 2 24cm Cat 4 SL	ROG2-24-4-SL
759324420	RUSH ROGUE 2 24cm Cat 4 SR	ROG2-24-4-SR
759324500	RUSH ROGUE 2 24cm cat 5	ROG2-24-5-RU
759324510	RUSH ROGUE 2 24cm Cat 5 SL	ROG2-24-5-SL
759324520	RUSH ROGUE 2 24cm Cat 5 SR	ROG2-24-5-SR
759324600	RUSH ROGUE 2 24cm cat 6	ROG2-24-6-RU
759324610	RUSH ROGUE 2 24cm Cat 6 SL	ROG2-24-6-SL
759324620	RUSH ROGUE 2 24cm Cat 6 SR	ROG2-24-6-SR
759325200	RUSH ROGUE 2 25cm cat 2	ROG2-25-2-RU
759325210	RUSH ROGUE 2 25cm Cat 2 SL	ROG2-25-2-SL
759325220	RUSH ROGUE 2 25cm Cat 2 SR	ROG2-25-2-SR
759325300	RUSH ROGUE 2 25cm cat 3	ROG2-25-3-RU
759325310	RUSH ROGUE 2 25cm Cat 3 SL	ROG2-25-3-SL
759325320	RUSH ROGUE 2 25cm Cat 3 SR	ROG2-25-3-SR
759325400	RUSH ROGUE 2 25cm cat 4	ROG2-25-4-RU
759325410	RUSH ROGUE 2 25cm Cat 4 SL	ROG2-25-4-SL
759325420	RUSH ROGUE 2 25cm Cat 4 SR	ROG2-25-4-SR
759325500	RUSH ROGUE 2 25cm cat 5	ROG2-25-5-RU
759325510	RUSH ROGUE 2 25cm Cat 5 SL	ROG2-25-5-SL
759325520	RUSH ROGUE 2 25cm Cat 5 SR	ROG2-25-5-SR
759325600	RUSH ROGUE 2 25cm cat 6	ROG2-25-6-RU
759325610	RUSH ROGUE 2 25cm Cat 6 SL	ROG2-25-6-SL
759325620	RUSH ROGUE 2 25cm Cat 6 SR	ROG2-25-6-SR
759325700	RUSH ROGUE 2 25cm cat 7	ROG2-25-7-RU
759325710	RUSH ROGUE 2 25cm Cat 7 SL	ROG2-25-7-SL
759325720	RUSH ROGUE 2 25cm Cat 7 SR	ROG2-25-7-SR
759326200	RUSH ROGUE 2 26cm cat 2	ROG2-26-2-RU
759326210	RUSH ROGUE 2 26cm Cat 2 SL	ROG2-26-2-SL
759326220	RUSH ROGUE 2 26cm Cat 2 SR	ROG2-26-2-SR
759326300	RUSH ROGUE 2 26cm cat 3	ROG2-26-3-RU
759326310	RUSH ROGUE 2 26cm Cat 3 SL	ROG2-26-3-SL
759326320	RUSH ROGUE 2 26cm Cat 3 SR	ROG2-26-3-SR
759326400	RUSH ROGUE 2 26cm cat 4	ROG2-26-4-RU
759326410	RUSH ROGUE 2 26cm Cat 4 SL	ROG2-26-4-SL

759326420	RUSH ROGUE 2 26cm Cat 4 SR	ROG2-26-4-SR
759326500	RUSH ROGUE 2 26cm cat 5	ROG2-26-5-RU
759326510	RUSH ROGUE 2 26cm Cat 5 SL	ROG2-26-5-SL
759326520	RUSH ROGUE 2 26cm Cat 5 SR	ROG2-26-5-SR
759326600	RUSH ROGUE 2 26cm cat 6	ROG2-26-6-RU
759326610	RUSH ROGUE 2 26cm Cat 6 SL	ROG2-26-6-SL
759326620	RUSH ROGUE 2 26cm Cat 6 SR	ROG2-26-6-SR
759326700	RUSH ROGUE 2 26cm cat 7	ROG2-26-7-RU
759326710	RUSH ROGUE 2 26cm Cat 7 SL	ROG2-26-7-SL
759326720	RUSH ROGUE 2 26cm Cat 7 SR	ROG2-26-7-SR
759327300	RUSH ROGUE 2 27cm cat 3	ROG2-27-3-RU
759327310	RUSH ROGUE 2 27cm Cat 3 SL	ROG2-27-3-SL
759327320	RUSH ROGUE 2 27cm Cat 3 SR	ROG2-27-3-SR
759327400	RUSH ROGUE 2 27cm cat 4	ROG2-27-4-RU
759327410	RUSH ROGUE 2 27cm Cat 4 SL	ROG2-27-4-SL
759327420	RUSH ROGUE 2 27cm Cat 4 SR	ROG2-27-4-SR
759327500	RUSH ROGUE 2 27cm cat 5	ROG2-27-5-RU
759327510	RUSH ROGUE 2 27cm Cat 5 SL	ROG2-27-5-SL
759327520	RUSH ROGUE 2 27cm Cat 5 SR	ROG2-27-5-SR
759327600	RUSH ROGUE 2 27cm cat 6	ROG2-27-6-RU
759327610	RUSH ROGUE 2 27cm Cat 6 SL	ROG2-27-6-SL
759327620	RUSH ROGUE 2 27cm Cat 6 SR	ROG2-27-6-SR
759327700	RUSH ROGUE 2 27cm cat 7	ROG2-27-7-RU
759327710	RUSH ROGUE 2 27cm Cat 7 SL	ROG2-27-7-SL
759327720	RUSH ROGUE 2 27cm Cat 7 SR	ROG2-27-7-SR
759327800	RUSH ROGUE 2 27cm cat 8	ROG2-27-8-RU
759327810	RUSH ROGUE 2 27cm Cat 8 SL	ROG2-27-8-SL
759327820	RUSH ROGUE 2 27cm Cat 8 SR	ROG2-27-8-SR
759327900	RUSH ROGUE 2 27cm cat 9	ROG2-27-9-RU
759327910	RUSH ROGUE 2 27cm Cat 9 SL	ROG2-27-9-SL
759327920	RUSH ROGUE 2 27cm Cat 9 SR	ROG2-27-9-SR
759328300	RUSH ROGUE 2 28cm cat 3	ROG2-28-3-RU
759328310	RUSH ROGUE 2 28cm Cat 3 SL	ROG2-28-3-SL
759328320	RUSH ROGUE 2 28cm Cat 3 SR	ROG2-28-3-SR
759328400	RUSH ROGUE 2 28cm cat 4	ROG2-28-4-RU
759328410	RUSH ROGUE 2 28cm Cat 4 SL	ROG2-28-4-SL
759328420	RUSH ROGUE 2 28cm Cat 4 SR	ROG2-28-4-SR
759328500	RUSH ROGUE 2 28cm cat 5	ROG2-28-5-RU
759328510	RUSH ROGUE 2 28cm Cat 5 SL	ROG2-28-5-SL
759328520	RUSH ROGUE 2 28cm Cat 5 SR	ROG2-28-5-SR
759328600	RUSH ROGUE 2 28cm cat 6	ROG2-28-6-RU
759328610	RUSH ROGUE 2 28cm Cat 6 SL	ROG2-28-6-SL
759328620	RUSH ROGUE 2 28cm Cat 6 SR	ROG2-28-6-SR
759328700	RUSH ROGUE 2 28cm cat 7	ROG2-28-7-RU
759328710	RUSH ROGUE 2 28cm Cat 7 SL	ROG2-28-7-SL
759328720	RUSH ROGUE 2 28cm Cat 7 SR	ROG2-28-7-SR
759328800	RUSH ROGUE 2 28cm cat 8	ROG2-28-8-RU
759328810	RUSH ROGUE 2 28cm Cat 8 SL	ROG2-28-8-SL
759328820	RUSH ROGUE 2 28cm Cat 8 SR	ROG2-28-8-SR
759328900	RUSH ROGUE 2 28cm cat 9	ROG2-28-9-RU
759328910	RUSH ROGUE 2 28cm Cat 9 SL	ROG2-28-9-SL
759328920	RUSH ROGUE 2 28cm Cat 9 SR	ROG2-28-9-SR
759329300	RUSH ROGUE 2 29cm cat 3	ROG2-29-3-RU
759329310	RUSH ROGUE 2 29cm Cat 3 SL	ROG2-29-3-SL
759329320	RUSH ROGUE 2 29cm Cat 3 SR	ROG2-29-3-SR
759329400	RUSH ROGUE 2 29cm cat 4	ROG2-29-4-RU
759329410	RUSH ROGUE 2 29cm Cat 4 SL	ROG2-29-4-SL
759329420	RUSH ROGUE 2 29cm Cat 4 SR	ROG2-29-4-SR
759329500	RUSH ROGUE 2 29cm cat 5	ROG2-29-5-RU
759329510	RUSH ROGUE 2 29cm Cat 5 SL	ROG2-29-5-SL
759329520	RUSH ROGUE 2 29cm Cat 5 SR	ROG2-29-5-SR
759329600	RUSH ROGUE 2 29cm cat 6	ROG2-29-6-RU
759329610	RUSH ROGUE 2 29cm Cat 6 SL	ROG2-29-6-SL
759329620	RUSH ROGUE 2 29cm Cat 6 SR	ROG2-29-6-SR
759329700	RUSH ROGUE 2 29cm cat 7	ROG2-29-7-RU
759329710	RUSH ROGUE 2 29cm Cat 7 SL	ROG2-29-7-SL
759329720	RUSH ROGUE 2 29cm Cat 7 SR	ROG2-29-7-SR
759329800	RUSH ROGUE 2 29cm cat 8	ROG2-29-8-RU
759329810	RUSH ROGUE 2 29cm Cat 8 SL	ROG2-29-8-SL
759329820	RUSH ROGUE 2 29cm Cat 8 SR	ROG2-29-8-SR
759329900	RUSH ROGUE 2 29cm cat 9	ROG2-29-9-RU
759329910	RUSH ROGUE 2 29cm Cat 9 SL	ROG2-29-9-SL
759329920	RUSH ROGUE 2 29cm Cat 9 SR	ROG2-29-9-SR