

KINTERRA 3.0

Kategorivelger i forhold til aktivitetsnivå, brukervekt og størrelse									
Vekt / Aktivitetsnivå	44-52 kg	53-59 kg	60-68 kg	69-77 kg	78-88 kg	89-100 kg	101-116 kg	117-132 kg	133-150 kg
Lavt	1	1	2	3	4	5	6	7	8
Moderat	1	2	3	4	5	6	7	8	-

Art.nr.	Modell	Lev.art.nr.
721301122	Kinterra 3.0 Foot & Ankle System, Cat. 1 Size 22cm SL	RM3-00-1A22-SL
721301123	Kinterra 3.0 Foot & Ankle System, Cat. 1 Size 23cm SL	RM3-00-1A23-SL
721301124	Kinterra 3.0 Foot & Ankle System, Cat. 1 Size 24cm SL	RM3-00-1A24-SL
721301125	Kinterra 3.0 Foot & Ankle System, Cat. 1 Size 25cm SL	RM3-00-1A25-SL
721301126	Kinterra 3.0 Foot & Ankle System, Cat. 1 Size 26cm SL	RM3-00-1A26-SL
721301127	Kinterra 3.0 Foot & Ankle System, Cat. 1 Size 27cm SL	RM3-00-1A27-SL
721301128	Kinterra 3.0 Foot & Ankle System, Cat. 1 Size 28cm SL	RM3-00-1A28-SL
721301222	Kinterra 3.0 Foot & Ankle System, Cat. 1 Size 22cm SR	RM3-00-1A22-SR
721301223	Kinterra 3.0 Foot & Ankle System, Cat. 1 Size 23cm SR	RM3-00-1A23-SR
721301224	Kinterra 3.0 Foot & Ankle System, Cat. 1 Size 24cm SR	RM3-00-1A24-SR
721301225	Kinterra 3.0 Foot & Ankle System, Cat. 1 Size 25cm SR	RM3-00-1A25-SR
721301226	Kinterra 3.0 Foot & Ankle System, Cat. 1 Size 26cm SR	RM3-00-1A26-SR
721301227	Kinterra 3.0 Foot & Ankle System, Cat. 1 Size 27cm SR	RM3-00-1A27-SR
721301228	Kinterra 3.0 Foot & Ankle System, Cat. 1 Size 28cm SR	RM3-00-1A28-SR
721302122	Kinterra 3.0 Foot & Ankle System, Cat. 2 Size 22cm SL	RM3-00-2A22-SL
721302123	Kinterra 3.0 Foot & Ankle System, Cat. 2 Size 23cm SL	RM3-00-2A23-SL
721302124	Kinterra 3.0 Foot & Ankle System, Cat. 2 Size 24cm SL	RM3-00-2A24-SL
721302125	Kinterra 3.0 Foot & Ankle System, Cat. 2 Size 25cm SL	RM3-00-2A25-SL
721302126	Kinterra 3.0 Foot & Ankle System, Cat. 2 Size 26cm SL	RM3-00-2A26-SL
721302127	Kinterra 3.0 Foot & Ankle System, Cat. 2 Size 27cm SL	RM3-00-2A27-SL
721302128	Kinterra 3.0 Foot & Ankle System, Cat. 2 Size 28cm SL	RM3-00-2A28-SL
721302222	Kinterra 3.0 Foot & Ankle System, Cat. 2 Size 22cm SR	RM3-00-2A22-SR
721302223	Kinterra 3.0 Foot & Ankle System, Cat. 2 Size 23cm SR	RM3-00-2A23-SR
721302224	Kinterra 3.0 Foot & Ankle System, Cat. 2 Size 24cm SR	RM3-00-2A24-SR
721302225	Kinterra 3.0 Foot & Ankle System, Cat. 2 Size 25cm SR	RM3-00-2A25-SR
721302226	Kinterra 3.0 Foot & Ankle System, Cat. 2 Size 26cm SR	RM3-00-2A26-SR
721302227	Kinterra 3.0 Foot & Ankle System, Cat. 2 Size 27cm SR	RM3-00-2A27-SR
721302228	Kinterra 3.0 Foot & Ankle System, Cat. 2 Size 28cm SR	RM3-00-2A28-SR
721303122	Kinterra 3.0 Foot & Ankle System, Cat. 3 Size 22cm SL	RM3-00-3A22-SL
721303123	Kinterra 3.0 Foot & Ankle System, Cat. 3 Size 23cm SL	RM3-00-3A23-SL
721303124	Kinterra 3.0 Foot & Ankle System, Cat. 3 Size 24cm SL	RM3-00-3A24-SL
721303125	Kinterra 3.0 Foot & Ankle System, Cat. 3 Size 25cm SL	RM3-00-3A25-SL
721303126	Kinterra 3.0 Foot & Ankle System, Cat. 3 Size 26cm SL	RM3-00-3A26-SL
721303127	Kinterra 3.0 Foot & Ankle System, Cat. 3 Size 27cm SL	RM3-00-3A27-SL
721303128	Kinterra 3.0 Foot & Ankle System, Cat. 3 Size 28cm SL	RM3-00-3A28-SL
721303222	Kinterra 3.0 Foot & Ankle System, Cat. 3 Size 22cm SR	RM3-00-3A22-SR
721303223	Kinterra 3.0 Foot & Ankle System, Cat. 3 Size 23cm SR	RM3-00-3A23-SR
721303224	Kinterra 3.0 Foot & Ankle System, Cat. 3 Size 24cm SR	RM3-00-3A24-SR
721303225	Kinterra 3.0 Foot & Ankle System, Cat. 3 Size 25cm SR	RM3-00-3A25-SR
721303226	Kinterra 3.0 Foot & Ankle System, Cat. 3 Size 26cm SR	RM3-00-3A26-SR
721303227	Kinterra 3.0 Foot & Ankle System, Cat. 3 Size 27cm SR	RM3-00-3A27-SR
721303228	Kinterra 3.0 Foot & Ankle System, Cat. 3 Size 28cm SR	RM3-00-3A28-SR

721307222	Kinterra 3.0 Foot & Ankle System, Cat. 7 Size 22cm SR	RM3-00-7A22-SR
721307223	Kinterra 3.0 Foot & Ankle System, Cat. 7 Size 23cm SR	RM3-00-7A23-SR
721307224	Kinterra 3.0 Foot & Ankle System, Cat. 7 Size 24cm SR	RM3-00-7A24-SR
721307225	Kinterra 3.0 Foot & Ankle System, Cat. 7 Size 25cm SR	RM3-00-7A25-SR
721307226	Kinterra 3.0 Foot & Ankle System, Cat. 7 Size 26cm SR	RM3-00-7A26-SR
721307227	Kinterra 3.0 Foot & Ankle System, Cat. 7 Size 27cm SR	RM3-00-7A27-SR
721307228	Kinterra 3.0 Foot & Ankle System, Cat. 7 Size 28cm SR	RM3-00-7A28-SR
721308122	Kinterra 3.0 Foot & Ankle System, Cat. 8 Size 22cm SL	RM3-00-8A22-SL
721308123	Kinterra 3.0 Foot & Ankle System, Cat. 8 Size 23cm SL	RM3-00-8A23-SL
721308124	Kinterra 3.0 Foot & Ankle System, Cat. 8 Size 24cm SL	RM3-00-8A24-SL
721308125	Kinterra 3.0 Foot & Ankle System, Cat. 8 Size 25cm SL	RM3-00-8A25-SL
721308126	Kinterra 3.0 Foot & Ankle System, Cat. 8 Size 26cm SL	RM3-00-8A26-SL
721308127	Kinterra 3.0 Foot & Ankle System, Cat. 8 Size 27cm SL	RM3-00-8A27-SL
721308128	Kinterra 3.0 Foot & Ankle System, Cat. 8 Size 28cm SL	RM3-00-8A28-SL
721308222	Kinterra 3.0 Foot & Ankle System, Cat. 8 Size 22cm SR	RM3-00-8A22-SR
721308223	Kinterra 3.0 Foot & Ankle System, Cat. 8 Size 23cm SR	RM3-00-8A23-SR
721308224	Kinterra 3.0 Foot & Ankle System, Cat. 8 Size 24cm SR	RM3-00-8A24-SR
721308225	Kinterra 3.0 Foot & Ankle System, Cat. 8 Size 25cm SR	RM3-00-8A25-SR
721308226	Kinterra 3.0 Foot & Ankle System, Cat. 8 Size 26cm SR	RM3-00-8A26-SR
721308227	Kinterra 3.0 Foot & Ankle System, Cat. 8 Size 27cm SR	RM3-00-8A27-SR
721308228	Kinterra 3.0 Foot & Ankle System, Cat. 8 Size 28cm SR	RM3-00-8A28-SR