

EVERYDAY

Freedom Agilix™

Maximum shock value

- Vertical shock absorption and carbon fiber flexibility optimized for low to moderate impact K3 ambulators; reduces socket shear forces and improves overall comfort
- Multi-axial function increases ground compliance and stability; streamlined, split-toe design results in improved long-term durability, saving time and money
- Full length dynamic heel attached at the toe ensures a seamless roll-over, for superior comfort and a more symmetrical gait
- Lower profile and lighter weight than the competitive alternative, allowing users to do more and feel less tired at the end of the day
- Rated to 166 kg
- 36-month warranty



Freedom Agilix (F15)



Vertical Shock




Split Toe



Sandal Toe

Freedom Agilix™



Note: Build height measurements represent a fully assembled unit.

Product ID: F15

Sizes: 22-31 cm

Sandal toe available for sizes 22-28 cm*

User Weight Rating: 166 kg

Heel Height: 10 mm

Connector: Male pyramid

Warranty: Foot module 36 months,
Cosmetic foot shell 6 months

*When ordering the Sandal Toe option, please contact Customer Service for lead time.

Accessories

- Foot Shell Removal Tool: (ACC-00-103)
- Spectra Sock: (ACC-00-101xx)
30 cm, 35 cm, 40 cm
- Stiffening Bumpers:
 - 22-25 cm (KIT-00-11363)
 - 26-28 cm (KIT-00-11364)
 - 29-31 cm (KIT-00-11365)

Note: Spectra socks and heel bumpers are included in the delivery.

Build Height & Weight

Size (cm)	Height (mm)	Weight (g)
22-25	143	553
26-28	161	609
29-31	172	772

Product weight includes foot shell. Measurements are approximate and will vary due to category selection and manufacturing variability.

Stiffness Category Selection

Impact Level		Weight (kg)
Low	Moderate	
1	1	44-52
1	2	53-59
2	3	60-68
3	4	69-77
4	5	78-88
5	6	89-100
6	7	101-116
7	8	117-130
8	9	131-147
9		148-166

Complete Prosthetic Foot (incl. foot shell) Part Number

F15-00-A--

Toe: 0-regular, 1-sandal
 Category: 1-9
 Size: 22-31
 Version: M-mid profile, S-sandal toe
 Side: L-left, R-right
 Generation: 3-regular toe, 2-sandal toe
 Skin Tone: L-light, M-medium, D-dark

Foot Module Part Number

F15-00-0A-

Category: 1-9
 Size: 22-31
 Toe Configuration: RU-regular toe, universal; SR-sandal toe, right; SL-sandal toe, left

Foot Shell Part Number

FTC--4-

Generation: 2-3
 Profile: M-mid, H-high
 Cap: 0=yes, 1=no
 Size: 22-31
 Skin Tone: L-light, M-medium, D-dark
 Toe: R-regular, S-sandal
 Side: L-left, R-right

See the Cosmetic Foot Shell page for more detailed ordering instructions.