



1

2

3

4



SPORT

# Freestyle™ Swim

## It's push button easy!

- Users can now walk proudly and comfortably to and from the swim area
- Adjust to 70° plantar flexion for swimming with the simple push of a button
- Utilize a flipper to gain swimming speed and decrease energy expenditure
- Balance the body and enhance swimming efficiency
- Rated to 116 kg
- 12-month warranty



Freestyle Swim (LP2-W2)




EnduraCore®  
Technology

# Freestyle™ Swim



SPORT



Note: Build height measurements represent a fully assembled unit.

**Product ID:** LP2-W2

**Sizes:** 22-31 cm

**User Weight Rating:** 116 kg

**Heel Height:** 10 mm

**Connector:** Male pyramid

**Warranty:** Foot module 36 months,  
Ankle 12 months, Cosmetic foot shell 6 months  
Freestyle Swim can be delivered with a special  
shaped cover, please contact Customer Service.

## Accessories

- Protective Cover for Ankle: (KIT-00-200)
- Foot Shell Removal Tool: (ACC-00-103)
- Spectra Sock: (ACC-00-101xx)  
30 cm, 35 cm, 40 cm
- Stiffening Bumpers:
  - 22-25cm (KIT-00-11100-00)
  - 26-28cm (KIT-00-11300-00)
  - 29-31cm (KIT-00-11500-00)

Note: Spectra socks and heel bumpers are included in the delivery.

## Build Height & Weight

Size (cm)	Height (mm)	Weight (g)
22-25	102	680
26-28	107	790
29-31	117	945

Product weight includes foot shell. Measurements are approximate and will vary due to category selection and manufacturing variability.

## Stiffness Category Selection

Impact Level			Weight (kg)
Low	Moderate	High	
1	1	2	44-52
1	2	3	53-59
2	3	4	60-68
3	4	5	69-77
4	5	6	78-88
5	6		89-100
6			101-116

## Complete Prosthetic Foot (incl. foot shell) Part Number

LP2-W2-1  A  -M  -3

Category 1-6      Size 22-31      Side L-left R-right      Skin Tone L-light M-medium D-dark

### Foot Module Part Number

LP2-W2-0  A  -RU

Category 1-6      Size 22-31

### Foot Shell Part Number

FTC-3M-1   4-R

Size 22-31      Skin Tone L-light M-medium D-dark      Side L-left R-right

See the Cosmetic Foot Shell page for more detailed ordering instructions.