



1

2

3

4



SPECIALTY

# LP Symes<sup>®</sup>

## Designed for the demanding needs of a Symes amputee

- Engineered to withstand the high forces of a Symes amputee, while remaining flexible enough to offer superior comfort
- Remarkably lightweight—averaging 295 g
- Direct lamination cup keeps build height low
- Rated up to 166 kg
- 36-month warranty




LP Symes (LP2)



EnduraCore<sup>®</sup>  
Technology

# LP Symes®



Note: Build height measurements represent a fully assembled unit.

**Product ID:** LP2

**Sizes:** 22-31 cm

**User Weight Rating:** 166 kg

**Heel Height:** 10 mm

**Connector:** Direct lamination cup

**Warranty:** Foot module 36 months,  
Cosmetic foot shell 6 months

## Accessories

- Direct Lamination Cup: (ACC-00-162)
- Foot Shell Removal Tool: (ACC-00-103)
- Spectra Sock: (ACC-00-101xx)  
30 cm, 35 cm, 40 cm
- Stiffening Bumpers:
  - 22-25 cm (KIT-00-111)
  - 26-28 cm (KIT-00-113)
  - 29-31 cm (KIT-00-115)

Note: Spectra socks and heel bumpers are included in the delivery.

## Build Height & Weight

Size (cm)	Height (mm)	Weight (g)
22-25	43	420
26-28	48	545
29-31	53	670

Product weight includes foot shell. Measurements are approximate and will vary due to category selection and manufacturing variability.

## Stiffness Category Selection

Impact Level			Weight (kg)
Low	Moderate	High	
1	1	2	44-52
1	2	3	53-59
2	3	4	60-68
3	4	5	69-77
4	5	6	78-88
5	6	7	89-100
6	7	8	101-116
7	8	9	117-130
8	9		131-147
9			148-166

## Complete Prosthetic Foot (incl. foot shell) Part Number

LP2-00-1  A  -M  -3

Category 1-9      Size 22-31      Side L-left R-right      Skin Tone L-light M-medium D-dark

### Foot Module Part Number

LP2-00-0  A  -RU

Category 1-9      Size 22-31

### Foot Shell Part Number

FTC-3M-1   4-R

Size 22-31      Skin Tone L-light M-medium D-dark      Side L-left R-right

See the Cosmetic Foot Shell page for more detailed ordering instructions.