

SPECIALTY

Pacifica[®] & Pacifica[®] LP

The perfect balance of flexibility and durability

- Choice of two low profile designs ideal for amputees with long residual limbs or those needing more build space
- LP design requires less than 50 mm of clearance
- Lightweight and durable
- Foot lies within foot shell so cosmetic finishing is simplified
- Rated up to 166 kg
- 36-month warranty



Pacifica (FS2)

Pacifica LP (FS4)



Split Toe




Sandal Toe



EnduraCore[®]
Technology

Pacifica® & Pacifica® LP



Note: Build height measurements represent a fully assembled unit.

Product ID: Pacifica (FS2); Pacifica LP (FS4)

Sizes: Pacifica 24-31 cm; Pacifica LP 22-31 cm

Sandal toe available for sizes 22-28 cm (FS2 24-28 cm)*

User Weight Rating: 166 kg

Heel Height: 10 mm

Connector: Male pyramid

Warranty: Foot module 36 months,
Cosmetic foot shell 6 months

*When ordering the Sandal Toe option, please contact Customer Service for lead time.

Accessories

- Foot Shell Removal Tool: (ACC-00-103)
- Spectra Sock: (ACC-00-101xx)
30 cm, 35 cm, 40 cm
- Stiffening Bumpers:

Pacifica (FS2):

- 24-25 cm (KIT-00-11343)
- 26-28 cm (KIT-00-11344)
- 29-31 cm (KIT-00-11345)

Pacifica LP (FS4):

- 22-25 cm (KIT-00-11323)
- 26-28 cm (KIT-00-11324)
- 29-31 cm (KIT-00-11325)

Note: Spectra socks and heel bumpers are included in the delivery.

Build Height & Weight

| | Size (cm) | Height (mm) | Weight (g) |
|-------------|-----------|-------------|------------|
| Pacifica | 24-25 | 59 | 485 |
| | 26-28 | 61 | 610 |
| | 29-31 | 63 | 720 |
| Pacifica LP | 22-25 | 47 | 465 |
| | 26-28 | 50 | 570 |
| | 29-31 | 60 | 760 |

Product weight includes foot shell. Measurements are approximate and will vary due to category selection and manufacturing variability.

Stiffness Category Selection

| Impact Level | | | Weight (kg) |
|--------------|----------|------|-------------|
| Low | Moderate | High | |
| 1 | 1 | 2 | 44-52 |
| 1 | 2 | 3 | 53-59 |
| 2 | 3 | 4 | 60-68 |
| 3 | 4 | 5 | 69-77 |
| 4 | 5 | 6 | 78-88 |
| 5 | 6 | 7 | 89-100 |
| 6 | 7 | 8 | 101-116 |
| 7 | 8 | 9 | 117-130 |
| 8 | 9 | | 131-147 |
| 9 | | | 148-166 |

Complete Prosthetic Foot (incl. foot shell) Part Number

_ _ _ -00- _ _ A _ - _ _ _ _

| | | | | | | | |
|------------|-----------|----------|-------|------------------|---------|---------------|-----------|
| Product ID | Toe | Category | Size | Version | Side | Generation | Skin Tone |
| FS2 | 0-regular | 1-9 | 24-31 | M-mid | L-left | 3-regular toe | L-light |
| FS4 | 1-sandal | Pacifica | 22-31 | Pacifica profile | R-right | 2-sandal toe | M-medium |
| | | | | S-sandal | | | D-dark |
| | | | | Pacifica LP | toe | | |

Foot Module Part Number

_ _ _ -00-0 _ A _ - _

| | | | |
|-------------------|----------|-------------|----------------------|
| Product ID | Category | Size | Toe Configuration |
| FS2 (Pacifica) | 1-9 | 24-31 (FS2) | RU-regular toe, |
| FS4 (Pacifica LP) | | 22-31 (FS4) | universal |
| | | | SR-sandal toe, right |
| | | | SL-sandal toe, left |

Foot Shell Part Number

FTC - _ _ - _ _ _ _ 4 - _ _

| | | | | | | |
|------------|-----------------|---------------|-------|-------------------------------|-----------------------|-------------------|
| Generation | Profile | Cap | Size | Skin Tone | Toe | Side |
| 2-3 | M-mid H-high | 0-yes 1-no | 22-31 | L-light M-medium D-dark | R-regular S-sandal | L-left R-right |

See the Cosmetic Foot Shell page for more detailed ordering instructions.