

Discover Meracus[®]

Pure composite
Pure motion



Pure composite. Pure motion.

Offer your patients the opportunity to move smoothly, enjoying a balanced energy return with every step they take.

Brought to you by the dedicated team behind Allard AFO, Meracus is a one-of-a-kind prosthetic foot crafted from 100% composite materials. Meracus is free from metal joints or any other extraneous components, making it a pure composite prosthetic foot. This design, in combination with the rocker sole, facilitates a smooth transition from heel strike to toe-off, ensuring a fluid, natural, and confident stride.

Let your patients discover a seamless gait experience with Meracus, and join us in our mission to provide Support for Better Life!



"Patients love the light weight of the foot and how compliant it is - even on uneven terrain."

- Prosthetist, Stockholm Sweden



100% composite precision for a seamless gait experience

Hybrid composite technology

A blend of fibers developed to maximize function throughout the gait cycle.

Enhanced roller shape

Heel and forefoot contours that improve loading response and stability.

Composite rivet technology

Heel and keel are joined with a unique composite technology that eliminates metal fasteners for improved flexibility and strength.

100% composite design

From the attachment to the toe for smooth transition and improved dynamic response.

Low profile J-shape

Longer composite curve for better shock absorption and energy return.


Extended heel lever design

Improved shock absorption and transition to mid-stance.

Meracus is an energy-storing, dynamic response prosthetic foot available in size 23-28, P3-P7.

"So soft and flexible, feels like my own foot!"

- 53 year old lower leg amputee





Support for Better Life

Everyone should be able to live their life to the fullest, regardless of their mobility challenges. With innovative solutions developed in close collaboration with healthcare professionals and patients, we strive to provide Support for Better Life.

allard

Allard International

Karbingatan 38, 254 67 Helsingborg, Sweden | +46 42 25 27 00 | info@allardint.com | allardint.com