

3204 XP Release Valve



The 3204 XP Release valve is an auto-expulsion with push button suction release in the same compact, easy to install component, and designed to be used with all suction suspension systems.

- Transtibial or Transfemoral users.
- Kit includes a flexible inner socket housing adapter ring.
- Threaded (M10) body to securely attach through the socket from the outside near the distal end.
- Auto expels air under positive pressure when entering a socket.
- Remove prosthesis with the push button release to allow air back into the socket.
- Can be retrofitted to an existing socket.
- Can be disassembled by hand.

VALVE INFORMATION

Part Number:	3204 (XP Release Value)
Expulsion Method:	Auto Expulsion
Amputation Level:	Transtibial, Transfemoral
Fabrication Method:	Lamination, Thermoplastic
Material:	Plastic
Product Dimension:	7/8" diameter
Weight:	10g
Warranty:	6 months from invoice date