

APC Estimate Build Height / Product Weight Chart

Part Number	Product Weight Rating	Description	Estimate Build Height (Inches)	Estimate Build Height (mm)	Weight (Including Screws) oz	Weight (Including Screws) grams
TITANIUM COMPONENTS						
TI-00MS	365 lbs / 166 kg	Titanium Rotatable Pyramid w/ Male Thread and Center Hole	0.24	6	1.7	48.19
TI-01MS	365 lbs / 166 kg	Titanium Rotatable Pyramid w/ Female Thread & Clamp	0.69	17.5	2.7	76.54
TI-100C	365 lbs / 166 kg	Titanium 30mm Tube Clamp Adapter	1.81	46	3	85.05
TI-100C-LONG	365 lbs / 166 kg	Titanium 30mm Tube Clamp Adapter LONG (84mm)	3.27	83	3.8	107.73
TI-100C-SH	365 lbs / 166 kg	Titanium 30mm Tube Clamp Adapter SHORT w/ Inside KNURL	1.5	38	2.8	79.38
TI-10TC	365 lbs / 166 kg	10" 30mm Pylon w/ Titanium Receiver Adapter	10.48	266.21	7	198.45
TI-10TC-I	365 lbs / 166 kg	10" 30mm Pylon w/ Internal Titanium Receiver Adapter	9.5	241.32	7.1	201.28
TI-105C	365 lbs / 166 kg	Titanium 30mm Tube Clamp Adapter w/ 5mm Offset Receiver	1.88	47.75	3.8	107.73
TI-110AT	365 lbs / 166 kg	Titanium 30mm Pylon (LONG) w/ Receiver Adapter - Overall Length = 3.8"	0.65	16.51	4	107.73
TI-110SH	365 lbs / 166 kg	Titanium 30mm Pylon (SHORT) w/ Receiver Adapter - Overall Length of 2.06" Accommodates 1 - TI-100C	0.75	19	2.7	76.54
TI-112-P	365 lbs / 166 kg	Titanium 30mm Pylon (SHORT) w/ Pyramid Adapter - Overall Length of 2.46" Accommodates 1 - TI-100C	0.63	16	3.1	87.88
TI-115DM	365 lbs / 166 kg	Titanium Dual Pyramid Adapter	0.6	15.24	4	113.4
TI-120P	365 lbs / 166 kg	Titanium 30mm Tube Clamp w/ Pyramid Adapter	1.26	32	3.1	87.88
TI-125P	365 lbs / 166 kg	Titanium Pyramid / Receiver Adapter	0.7	17.78	2.6	73.71
TI-132D	365 lbs / 166 kg	Titanium 32mm Dual Receiver Adapter	1.7	43.18	3.5	99.22
TI-145D	365 lbs / 166 kg	Titanium 45mm Dual Receiver Adapter	1.75	44.45	4	113.4
TI-160D	365 lbs / 166 kg	Titanium 60mm Dual Receiver Adapter	2.25	57.15	4.5	127.57
TI-175D	365 lbs / 166 kg	Titanium 75mm Dual Receiver Adapter	2.85	72.39	5.5	155.92
TI-17TC	365 lbs / 166 kg	17" 30mm Pylon w/ Titanium Receiver Adapter	17.48	444	10	283.49
TI-17TC-I	365 lbs / 166 kg	17" 30mm Pylon w/ Internal Titanium Receiver Adapter	16.5	419.12	10.9	309.01
TI-205C	365 lbs / 166 kg	Titanium 5 Degree Angled Tube Clamp Adapter	1.93	49	3	85.05
TI-210C	365 lbs / 166 kg	Titanium 10 Degree Angled Tube Clamp Adapter	2.05	52	3.2	90.72
TI-215C	365 lbs / 166 kg	Titanium 15 Degree Angled Tube Clamp Adapter	2.13	54	3.4	96.39
TI-220C	365 lbs / 166 kg	Titanium 20 Degree Angled Tube Clamp Adapter	2.24	57	3.6	102.06
TI-400A	365 lbs / 166 kg	Titanium Offset 4 Hole Pyramid Adapter (Axial)	0.31	8	2.9	82.21
TI-400B	365 lbs / 166 kg	Titanium 4 Hole Pyramid Adapter	0.31	8	1.75	49.61
TI-400BIL	365 lbs / 166 kg	Titanium 4 Hole Pyramid Adapter (ROUND Base)	0.31	8	2.1	59.53
TI-400D	365 lbs / 166 kg	Titanium Offset 4 Hole Pyramid Adapter (Diagonal)	0.31	8	2.9	82.21

Part Number	Product Weight Rating	Description	Estimate Build Height (Inches)	Estimate Build Height (mm)	Weight (Including Screws) oz	Weight (Including Screws) grams
TI-400RM	365 lbs / 166 kg	Titanium 4 Hole Rotatable Pyramid Adapter	0.38	9.65	2.1	59.53
TI-406PAT	365 lbs / 166 kg	Titanium Rotatable Pyramid Adapter w/ Threaded 4 Hole Base	0.38	9.65	2.5	7.87
TI-420L	365 lbs / 166 kg	Titanium Round 4 Hole Threaded Plate (.188" Thick)	0.16	4	1.1	31.18
TI-425BH	365 lbs / 166 kg	Titanium 4 Hole Pyramid Adapter w/ Center Hole	0.31	8	1.75	49.61
TI-450BF	365 lbs / 166 kg	Titanium 4 Hole Receiver Adapter	0.7	17.78	3	85.05
TI-455RF	365 lbs / 166 kg	Titanium 4 Hole Rotatable Receiver Adapter	0.55	13.97	2.5	70.87
TI-456PAT	365 lbs / 166 kg	Titanium Rotatable Receiver Adapter w/ 4 Hole Threaded Base	0.55	13.97	2.5	70.87
TI-460B	365 lbs / 166 kg	Titanium 4 Hole Clamping Adapter w/ Center Thread (BASE ONLY)	0.38	9.65	2.2	62.37
TI-460TB	365 lbs / 166 kg	Titanium 4 Hole Clamping Adapter w/ Center Thread (BASE ONLY) All 4 Holes Counter Sunk	0.31	8	1.9	53.86
TI-460TRM	365 lbs / 166 kg	Titanium 4 Hole Clamping Adapter (Center Thread) w/ Rotatable Pyramid - All 4 Holes Counter Sunk	0.55	14	3.5	99.22
TI-465RF	365 lbs / 166 kg	Titanium 4 Hole Clamping Adapter (Center Thread) w/ Rotatable Receiver	0.83	21	3.9	110.56
TI-465TRF	365 lbs / 166 kg	Titanium 4 Hole Clamping Adapter (Center Thread) w/ Rotatable Receiver - All 4 Holes Counter Sunk	0.83	21	3.8	107.73
TI-4PRWH	365 lbs / 166 kg	Titanium 4 Hole Pyramid with Center Hole (Round Base w/ Rotation Slots)	0.27	7	2.5	70.87
TI-4RR	365 lbs / 166 kg	Titanium 4 Hole Receiver (Round Base w/ Rotation Slots)	0.63	16	3.4	96.38
TI-500P	365 lbs / 166 kg	Titanium Foot Pyramid Adapter (Bolt and Pins Included)	0.25	6.35	4	113.4
TI-50FS	365 lbs / 166 kg	Titanium Rotatable Receiver w/ Male Thread	0.51	13	1.9	53.86
TI-600AD	365 lbs / 166 kg	Titanium 4 Hole Alignable Pyramid Adapter	0.78	19.81	6.3	178.6
TI-600-SLA-PM	365 lbs / 166 kg	Titanium (Sleeve System) Pyramid w/ Male Dovetail Slide Adapter	0.7	17.78	3.1	87.88
TI-60FS	365 lbs / 166 kg	Titanium Receiver w/ Female Thread	0.83	21.08	2.9	82.21
TI-610-SLA-TF	365 lbs / 166 kg	Titanium (Sleeve System) Male Thread w/ Female Dovetail Slide Adapter	0.64	16.26	2.8	79.38
TI-615-SLA-CF	365 lbs / 166 kg	Titanium (Sleeve System) 30mm Tube Clamp w/ Female Dovetail Slide Adapter	0.72	18.29	3.6	102.06
TI-620-SLA-FF	365 lbs / 166 kg	Titanium (Sleeve System) Dual Female Dovetail Slide Adapter	1.25	31.75	3.4	96.39
TI-630-SLA-RF	365 lbs / 166 kg	Titanium (Sleeve System) Receiver w/ Female Dovetail Slide Adapter	1.25	31.75	3.6	102.06
TI-640-SLA-PF	365 lbs / 166 kg	Titanium (Sleeve System) Pyramid w/ Female Dovetail Slide Adapter	0.8	20.32	3.3	93.55
TI-650-SLA-RM	365 lbs / 166 kg	Titanium (Sleeve System) Receiver w/ Male Dovetail Slide Adapter	0.68	17.27	3.6	102.06
TI-655-SLA-CM	365 lbs / 166 kg	Titanium (Sleeve System) 30mm Tube Clamp w/ Male Dovetail Slide Adapter	0.3	7.62	3.6	102.06
TI-675AR	365 lbs / 166 kg	Titanium Rotatable and Alignable Dual Receiver Adapter (7.5mm Adjustment)	1.4	35.56	6.7	189.94
TI-685AR	365 lbs / 166 kg	Titanium 4 Hole Alignable & Rotatable Receiver Adapter	1.38	35.05	7.3	206.95
TI-700LM	365 lbs / 166 kg	Stainless Steel 3 Prong Laminating Base (MALE THREAD) w/ Titanium Rotatable Pyramid Adapter	0.75	19.05	7.1	201.28
TI-740-B	365 lbs / 166 kg	Titanium 4 Prong Laminating Base (MALE THREAD)	0.47	12	2.1	59.53

Part Number	Product Weight Rating	Description	Estimate Build Height (Inches)	Estimate Build Height (mm)	Weight (Including Screws) oz	Weight (Including Screws) grams
TI-740L	365 lbs / 166 kg	Titanium 4 Prong Laminating Base w/ Fixed Pyramid Adapter	0.53	13.5	2	56.7
TI-740L-S	365 lbs / 166 kg	Titanium 4 Prong Laminating Adapter (SHORT)	0.51	13	2.9	82.21
TI-740-RM	365 lbs / 166 kg	Titanium 4 Prong Laminating Base (MALE THREAD) w/ Titanium Rotatable Pyramid Adapter	0.8	20.32	4.5	127.57
TI-745-RF	365 lbs / 166 kg	Titanium 4 Prong Laminating Base (MALE THREAD) w/ Titanium Rotatable Receiver Adapter	1.2	30.48	4.7	133.24
TI-750LM	365 lbs / 166 kg	Stainless Steel 3 Prong Laminating Base (MALE THREAD) w/ Titanium Rotatable Receiver Adapter	0.98	24.89	7	198.45
TI-800LM	365 lbs / 166 kg	Stainless Steel 3 Prong Laminating Base (Female Thread) w/ Titanium Rotatable Pyramid Adapter (Includes Lamination Cap)	0.55	14	5.3	150.25
TI-840LM	365 lbs / 166 kg	Stainless Steel 4 Prong Laminating Base (Female Thread) w/ Titanium Rotatable Pyramid Adapter (Includes Lamination Cap)	0.79	20	5.1	144.58
TI-845LM	365 lbs / 166 kg	Stainless Steel 4 Prong Laminating Base (Female Thread) w/ Titanium Rotatable Receiver Adapter (Includes Lamination Cap)	1.06	27	5.3	150.25
TI-850LM	365 lbs / 166 kg	Stainless Steel 3 Prong Laminating Base (Female Thread) w/ Titanium Rotatable Receiver Adapter (Includes Lamination Cap)	0.79	20	5.5	155.92
TI-CPF-400	365 lbs / 166 kg	Titanium 4 Hole Threaded Laminating Ring (Includes Lamination Cap)	0	0	3.3	93.55
TI-CPF-800	365 lbs / 166 kg	Titanium Laminating Pyramid Adapter (Includes Lamination Cap)	0.63	16	2.4	68.04
TI-TF900SET	365 lbs / 166 kg	Titanium Thermoplastic Socket Attachment Plate Set	0.49	12.45	4	113.4
STAINLESS STEEL COMPONENTS						
SS-00MS	365 lbs / 166 kg	Stainless Steel Rotatable Pyramid w/ Male Thread and Center Hole	0.24	6	3.3	93.55
SS-01MS	365 lbs / 166 kg	Stainless Steel Rotatable Pyramid w/ Female Thread & Clamp	0.6	15.24	4.8	136.08
SS-100C	365 lbs / 166 kg	Stainless Steel 30mm Tube Clamp Adapter	1.75	44.45	5	141.75
SS-10TC	365 lbs / 166 kg	10" 30mm Pylon w/ Stainless Steel Receiver Adapter	10.55	268	6.9	195.61
SS-132D	365 lbs / 166 kg	Stainless Steel 32mm Dual Receiver Adapter	1.7	43.18	4.9	138.91
SS-145D	365 lbs / 166 kg	Stainless Steel 45mm Dual Receiver Adapter	1.75	44.45	6	170.1
SS-160D	365 lbs / 166 kg	Stainless Steel 60mm Dual Receiver Adapter	2.25	57.15	7.4	209.79
SS-175D	365 lbs / 166 kg	Stainless Steel 75mm Dual Receiver Adapter	2.85	72.39	8.7	246.64
SS-17TC	365 lbs / 166 kg	17" 30mm Pylon w/ Stainless Steel Receiver Adapter	17.56	446	10.2	289.17
SS-300LM	365 lbs / 166 kg	Stainless Steel 3 Prong Laminating Base (MALE THREAD)	0.28	7	4.4	124.74
SS-400B	365 lbs / 166 kg	Stainless Steel 4 Hole Pyramid Adapter	0.28	7.11	3.3	93.55
SS-425BH	365 lbs / 166 kg	Stainless Steel 4 Hole Pyramid Adapter w/Center Hole	0.28	7.11	3.2	90.72
SS-460TB	365 lbs / 166 kg	Stainless Steel 4 Hole Clamping Adapter w/ Female Thread (BASE ONLY)	0.35	9	2.5	70.87
SS-500P	365 lbs / 166 kg	Stainless Steel Foot Pyramid Adapter (Bolt and Pins Included)	0.25	6.35	4.4	124.74
SS-50FS	365 lbs / 166 kg	Stainless Steel Rotatable Receiver w/ Male Thread	0.51	13	3	85.05
SS-60FS	365 lbs / 166 kg	Stainless Steel Rotatable Receiver w/ Female Thread & Clamp	0.91	23	4.2	119.07
SS-700LM	365 lbs / 166 kg	Stainless Steel 3 Prong Laminating Base (MALE THREAD) w/ Stainless Steel Rotatable Pyramid Adapter	0.75	19.05	9.25	262.23

Part Number	Product Weight Rating	Description	Estimate Build Height (Inches)	Estimate Build Height (mm)	Weight (Including Screws) oz	Weight (Including Screws) grams
SS-740L	365 lbs / 166 kg	Stainless Steel 4 Prong Laminating Base w/ Fixed Pyramid Adapter	0.25	6.35	3.5	99.22
SS-750LM	365 lbs / 166 kg	Stainless Steel 3 Prong Laminating Base (MALE THREAD) w/ Stainless Steel Rotatable Receiver Adapter	0.98	24.89	9	255.15
SS-800BASE	365 lbs / 166 kg	Stainless Steel 3 Prong Laminating Base w/ Female Thread (Includes Lamination Cap)	0.31	8	3.6	102.06
SS-800LM	365 lbs / 166 kg	Stainless Steel 3 Prong Laminating Base (Female Thread) w/ Stainless Steel Rotatable Pyramid Adapter (Includes Lamination Cap)	0.55	14	6.6	187.11
SS-810LM	365 lbs / 166 kg	Stainless Steel 3 Prong Laminating Base w/ Fixed Stainless Steel Pyramid Adapter (Includes Lamination Cap)	0.47	12	4.5	127.57
SS-840BASE	365 lbs / 166 kg	Stainless Steel 4 Prong Laminating Base w/ Female Thread (Includes Lamination Cap)	0.51	13	3.4	96.39
SS-840LM	365 lbs / 166 kg	Stainless Steel 4 Prong Laminating Base (Female Thread) w/ Stainless Steel Rotatable Pyramid Adapter (Includes Lamination Cap)	0.79	20	6.5	184.27
SS-845LM	365 lbs / 166 kg	Stainless Steel 4 Prong Laminating Base (Female Thread) w/ Stainless Steel Rotatable Receiver Adapter (Includes Lamination Cap)	1.06	27	6.4	181.44
SS-850LM	365 lbs / 166 kg	Stainless Steel 3 Prong Laminating Base (Female Thread) w/ Stainless Steel Rotatable Receiver Adapter (Includes Lamination Cap)	0.79	20	6.5	184.27
HEAVY DUTY COMPONENTS						
HD-100C	500 lbs / 228 kg	Heavy Duty Titanium 34mm Tube Clamp (Rated to 500 lbs)	2.05	52.07	5.4	153.09
HD-100C-A	500 lbs / 228 kg	Heavy Duty Titanium 34mm Tube Clamp w/30mm Receiver Adapter (Rated to 500 lbs)	0.6	15.24	4.3	121.9
HD-105C	500 lbs / 228 kg	Heavy Duty Titanium 34mm Tube Clamp w/5mm Offset (Rated to 500 lbs)	2	50.8	6	170.1
HD-110AT	500 lbs / 228 kg	Heavy Duty Titanium 34mm Pylon w/ Receiver Adapter (Rated to 500 lbs)	0.82	20.83	5.1	144.58
HD-110AT-LONG	500 lbs / 228 kg	Heavy Duty Titanium 34mm Pylon (LONG) w/ Receiver Adapter - LONG (PYLON = 2.00" - 3.21" Overall) (Rated to 500 lbs)	1.63	41.275	5.8	164.43
HD-132D	500 lbs / 228 kg	Heavy Duty Titanium 34mm Dual Receiver Adapter (Rated to 500 lbs)	1.75	44.45	7.4	209.79
HD-425BH	500 lbs / 228 kg	Heavy Duty Titanium 34mm Pyramid Adapter (Rated to 500 lbs)	0.5	12.7	4.6	130.41
HD-500P	500 lbs / 228 kg	Heavy Duty Titanium 34mm Foot Pyramid Adapter w/ Bolt (Rated to 500 lbs)	0.4	10.16	6	170.1
HD-740L	500 lbs / 228 kg	Heavy Duty Titanium 34mm 4 Prong Laminating Pyramid (Rated to 500 lbs)	0.65	16.51	7.2	204.12
HD-745L	500 lbs / 228 kg	Heavy Duty Titanium 34mm 4 Prong Laminating Receiver (Rated to 500 lbs)	0.75	19.05	6.6	187.11
HD-745L-A	500 lbs / 228 kg	Heavy Duty Titanium 4 Prong w/30mm Receiver (Rated to 500 lbs)	0.5	12.7	5.1	144.58
HD-AL-PYLON	500 lbs / 228 kg	Heavy Duty 34mm Aluminum Pylon 11.75" Length (Rated to 500 lbs)	11.75	298.45	14.8	419.57
ALUMINUM COMPONENTS						
AL-10T	365 lbs / 166 kg	10" 30mm Seamless Pylon	10	254	4.7	133.2427586
AL-17T	365 lbs / 166 kg	17" 30mm Seamless Pylon	17	431.8	8	226.796
AL-410-FLPL	365 lbs / 166 kg	Aluminum 10 Degree Flexion Plate	0.69	17.4625	6	170.1
AL-415-FLPL	365 lbs / 166 kg	Aluminum 15 Degree Flexion Plate	0.69	17.4625	6	170.1
AL-413-OS	365 lbs / 166 kg	Aluminum 4 Hole 13mm Offset Plate	0.25	6.35	3	85.05
AL-419-OS	365 lbs / 166 kg	Aluminum 4 Hole 19mm Offset Plate	0.25	6.35	3.4	96.39
AL-4X4-SLIDE	365 lbs / 166 kg	Aluminum 4 Hole 25.4mm Slide (Bolts Included)	0.5	12.7	4.5	127.57

Part Number	Product Weight Rating	Description	Estimate Build Height (Inches)	Estimate Build Height (mm)	Weight (Including Screws) oz	Weight (Including Screws) grams
AL-425-OS	365 lbs / 166 kg	Aluminum 4 Hole 25mm Offset Plate	0.25	6.35	3.5	99.22
AL-409-OS	365 lbs / 166 kg	Aluminum 4 Hole 9mm Offset Plate	0.25	6.35	2.9	82.21
AL-420L-SL	365 lbs / 166 kg	Aluminum 4 Hole Laminating Plate with Centered 6 Hole Shuttle Lock Relief	0	0	1.3	36.68
ALH-425BH	250 lbs / 114 kg	Aluminum 4 Hole Pyramid Adapter w/Center Hole (250 lb Rating)	0.34	8.64	2.1	59.53
AL-405-FLPL	365 lbs / 166 kg	Aluminum 5 Degree Flexion Plate	0.69	17.4625	6	170.1
AL-408-FLPL	365 lbs / 166 kg	Aluminum 8 Degree Flexion Plate	0.69	17.4625	6	170.1
ALH-500P	250 lbs / 114 kg	Aluminum Foot Pyramid Adapter (250 lb Rating)	0.34	8.64	4	113.4
AL-CPF-400	365 lbs / 166 kg	Aluminum Round 4 Hole Laminating Plate Anodized BLACK (5/8" Thick with 1/2" Tie Groove) Includes Screws	0	0	4.9	138.91
AL-420L-OS	365 lbs / 166 kg	Aluminum Round 4 Hole Laminating Plate w/offset Center & 6 Hole Shuttle Lock Relief	0	0	1.3	36.85
AL-420L-TH	365 lbs / 166 kg	Aluminum Round 4 Hole Laminating Plate with Threaded Center Hole (3/4 x 16)	0	0	2.5	70.87
AL-418B	365 lbs / 166 kg	Aluminum Round 4 Hole Spacer Plate (1/8")	0.13	3.18	0.5	14.17
AL-440L	365 lbs / 166 kg	Aluminum Round 4 Hole Threaded Laminating Ring (3/8" Thick with 1/8" Tie Groove)	0	0	2.4	68.04
AL-420L	365 lbs / 166 kg	Aluminum Round 4 Hole Threaded Lamination Plate	0	0	2.5	70.87
AL-431L	365 lbs / 166 kg	Aluminum Round 4 Hole Threaded Lamination Plate (.313)	0	0	3.8	107.73
AL-416B	365 lbs / 166 kg	Aluminum Square 4 Hole Spacer Plate (.08)	0.08	2.03	0.25	7.09
AL-420B	365 lbs / 166 kg	Aluminum Square 4 Hole Spacer Plate (.134)	0.19	4.7625	0.4	11.34
AL-430L	365 lbs / 166 kg	Aluminum Square 4 Hole Threaded Lamination Plate	0	0	2	56.7
AL-TF900SET	365 lbs / 166 kg	Aluminum Thermoplastic Socket Attachment Plate Set	0.49	12.45	4.2	119.07
PEDIATRIC/GERIATRIC COMPONENTS						
JR-100C	180 lbs / 81.65 kg	Pediatric / Geriatric Titanium 22mm Tube Clamp Adapter	1.55	39.37	2	56.7
JR-100C-A	180 lbs / 81.65 kg	Pediatric / Geriatric Titanium 22mm Tube Clamp Adapter w/30mm Receiver	1.75	44.45	3.4	96.39
JR-100C-LONG	180 lbs / 81.65 kg	Pediatric / Geriatric Titanium 22mm Tube Clamp Adapter LONG (2 5/8" Overall Height)	2.58	65.53	2	56.7
JR-100C-LONG-A	180 lbs / 81.65 kg	Pediatric / Geriatric Titanium 22mm Tube Clamp Adapter LONG w/30mm Adult Receiver (2 5/8" Overall Height)	2.63	66.675	2.8	79.38
JR-10T	180 lbs / 81.65 kg	Pediatric / Geriatric 10" 22mm Pylon	10	254	3.6	102.06
JR-10TC	180 lbs / 81.65 kg	Pediatric / Geriatric 10" 22mm Pylon w/Titanium Receiver	1.48	37.59	5.3	150.25
JR-10TC-A	180 lbs / 81.65 kg	Pediatric / Geriatric 10" 22mm Pylon w/30mm Adult Titanium Receiver	1.65	41.91	6.2	175.77
JR-10TC-I	180 lbs / 81.65 kg	Pediatric / Geriatric 10" 22mm Pylon w/Internal Titanium Receiver	10.25	260.35	4.7	133.24
JR-400A	180 lbs / 81.65 kg	Pediatric / Geriatric Titanium 4 Hole Offset Pyramid (AXIAL)	0.25	6.35	1.6	45.36
JR-400B	180 lbs / 81.65 kg	Pediatric / Geriatric Titanium 4 Hole Pyramid	0.25	6.35	1.3	36.85
JR-400D	180 lbs / 81.65 kg	Pediatric / Geriatric Titanium 4 Hole Offset Pyramid (DIAGONAL)	0.25	6.35	1.6	45.36

Part Number	Product Weight Rating	Description	Estimate Build Height (Inches)	Estimate Build Height (mm)	Weight (Including Screws) oz	Weight (Including Screws) grams
JR-416B	180 lbs / 81.65 kg	Pediatric / Geriatric Aluminum 4 Hole Spacer Plate (1/16" Thick)	0.06	1.52	0.01	0.28
JR-425BH	180 lbs / 81.65 kg	Pediatric / Geriatric Titanium 4 Hole Pyramid w/Center Hole	0.25	6.35	1.2	34.02
JR-450BF	180 lbs / 81.65 kg	Pediatric / Geriatric Titanium 4 Hole Receiver Adapter	0.7	17.78	1.9	53.86
JR-500P	180 lbs / 81.65 kg	Pediatric / Geriatric Foot Pyramid Adapter (Includes Bolt)	0.18	4.57	2	56.7
JR-AL-417B	180 lbs / 81.65 kg	Pediatric / Geriatric Aluminum 4 Hole Spacer (.134 Thick)	0.13	3.175	0.4	11.34
JR-AL-420L	180 lbs / 81.65 kg	Pediatric / Geriatric Aluminum Round 4 Hole Threaded Plate (M6 Screws Included)	0	0	2.4	68.04
JR-AL-420L-M5	180 lbs / 81.65 kg	Pediatric / Geriatric Aluminum Round 4 Hole Threaded Plate (M5 Screws Included)	0	0	2.2	62.37
JR-AL-431L-M5	180 lbs / 81.65 kg	Pediatric / Geriatric Round 4 Hole Threaded Plate .313 Thick (M5 Screws Included)	0	0	2.6	73.71
JR-AL-4X4-SLIDE	180 lbs / 81.65 kg	Pediatric / Geriatric Aluminum 4 Hole 18mm Slide (M5 Screws Included)	0.38	9.652	2.9	82.21