

Learn more about the L300 Foot Drop System  
by visiting our website [www.bioness.com](http://www.bioness.com).

JOIN US:



[www.youtube.com/user/BionessLiveOn](http://www.youtube.com/user/BionessLiveOn)



[www.twitter.com/BionessLiveOn](http://www.twitter.com/BionessLiveOn)



[www.facebook.com/BionessLiveOn](http://www.facebook.com/BionessLiveOn)



25103 Rye Canyon Loop

Valencia, CA 91355

800.211.9136

661.362.4850

[info@bioness.com](mailto:info@bioness.com)

[www.bioness.com](http://www.bioness.com)

1. Hausdorff JM, et al. *Am J Phys Med Rehab*. 2008;87(11):4-13.

2. Hausdorff JM, et al. *J Neurol Phys Ther*. 30(4):207.

Individual results vary. Consult with a qualified physician to determine if this product is right for you.

Contraindications, Adverse Reactions and Precautions are available on-line at [www.bioness.com](http://www.bioness.com) (also available in the L300 User's Guide).

NESS L300®, Intelli-Gait®, Intelli-Sense Gait Sensor®, Bioness®, the Bioness Logo and LiveOn® are trademarks of Bioness Inc. | [www.bioness.com](http://www.bioness.com)

**Rx Only**

© 2013 Bioness Inc.

Living with  
*Foot Drop?*



Regain  
*your*  
**Mobility**  
&  
Independence

NESS  
**L300®**



“A life changing experience for the **GOOD!**”

– Lorrie, L300 User

## Walk More **NATURALLY**

Foot Drop is a condition where the muscles in the foot are too weak to properly lift the foot and toes while walking. The award-winning L300® Foot Drop System helps to alleviate walking challenges that may result from:

- Stroke
- Traumatic brain injury
- Incomplete spinal cord injury
- Multiple sclerosis
- Cerebral Palsy

The L300 System is a functional electrical stimulation (FES) system that stimulates the nerves in the lower leg, activating muscles to lift the foot which allows you to walk more naturally and enjoy greater independence and a better quality of life. The system's gait sensing technology automatically detects your foot position, walking speed, and changes in terrain.

## Life Changing **BENEFITS**

The L300 System offers a variety of proven advantages to help people with foot drop:

- Provides a more natural movement when walking<sup>1</sup>
- Increase speed, stability and confidence<sup>1</sup>
- Reduce falls<sup>2</sup>
- Reeducate muscles
- Prevent muscle loss (often called atrophy)
- Maintain or increase range of motion in the ankle and foot
- Increase blood flow

Start your **JOURNEY**  
today with a  
**FREE SCREENING**

Call 1.800.211.9136, Option 2  
or visit our website at  
[www.bioness.com](http://www.bioness.com)

## **FREEDOM** with Wireless Design

The L300 System consists of three small components that communicate with each other wirelessly and adapt to changes in walking speed and terrain.

### **1** Ergonomic Leg Cuff

The comfortable leg cuff is free of control dials or switches to prevent accidentally changing your program settings and can be worn under most clothing.

### **2** Intelli-Sense Gait Sensor

Worn in the shoe, the Intelli-Sense Gait Sensor sends wireless signals to the leg cuff which stimulates the nerves to contract the muscles that lift your foot while walking.

### **3** Wireless Control Unit

A small hand held control unit communicates wirelessly with the system and allows you to safely and easily adjust the stimulation on the go.



## The Perfect **FIT**

The L300 System leg cuff is now available in two sizes to fit most adults and children:

Regular L300:

Fits calf circumference of 11 to 20 inches

Small L300:

Fits calf circumference of ~9 to 12 inches



**NESS**  
**L300®**